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Examine primary media plates daily for 7 days with dissecting microscope

Suspicious colonies are present.

No suspicious colonies present: **Culture negative**

Characteristics (on BCYE) can be:

- round, 1 – 2 mm, glistening colony at 3 days; larger with longer incubation
- cut-glass appearance, or with appearance of white pigment;
- green, or pink-purple, iridescence in young colonies, becoming opaque with longer incubation.
- autofluorescent under UV light demonstrated by some species other than *Legionella pneumophila* (Lp)

Subculture isolate to BCYE and BAP (presumptive test).

No growth on BAP; growth on BCYE:
Presumptive positive for *Legionella*

Growth on BAP; growth or no growth on BCYE: Not *Legionella*

Autofluorescence under UV light?

Yes – Orange, red, or blue-white autofluorescence:

No autofluorescence present but has typical Lp colonial morphology
Perform Latex Agglutination Test (LAT) : Lp serogroup 1 (Lp 1), Lp serogroups 2-14 (Lp2-14)

Orange or red:
Perform Direct
Fluorescent antibody (DFA):
L. rubrilucens, *L. erythra*,
Legionella sp. IFA

Blue-white:
Perform DFA;
L. dumoffii, *L. gormanii*,
L. bozemanii, and *L. anisa*

Positive Lp1 LAT
Perform Lp1 DFA

Negative Lp1 LAT
Positive Lp2-14 LAT
Perform DFA with mTech Lp1-14
And individual mTech reagents to Lp2 -6

Negative Lp1 LAT and Lp2-14 LAT
Perform *Legionella*.sp. FA:
Remel IFA or mTech Leg. spp.

If positive
with one:
Confirm ID as
L. rubrilucens, *L. erythra*,
or *Legionella* sp.,
not *pneumophila*.

If negative:
not *Legionella*.

If positive:
With one:
Confirm ID as
that specific
species.

If negative:
Perform
Legionella .sp FA
(Remel or mTech).
If (+): Confirm
ID as *Legionella* sp.,
not *pneumophila*.
If (-): then
not *Legionella*

If equivocal or cross-reacting
With more than one reagent:
Confirm as blue-white
fluorescing *Legionella* sp., not
pneumophila

If Positive:
Confirm
as Lp1.

If Negative:
Retest LAT.
Retest DFA
using
monoclonal
subtypes of
Lp1

If Positive with:
Lp2,Lp3,
Lp4 or Lp6,
and Lp1-14

Confirm as
that specific
serogroup.

Lp1-14 and
Lp5

Retest with
DFA Lp5
and Lp8 to
differentiate
and confirm.

Lp1-14
Only

Confirm
as *L. pneumophila*,
not serogroup 1-6

If Negative:
Retest LAT.

If FA Positive:
Confirm as
Legionella sp.,
not *pneumophila*

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